

Task Number #210097

Testing ordered > 19 JUL '26

Sample received > 23 JUL '26

Clientdrakonpeptides.com

SampleTesamorelin 10mg

ManufacturerUnknown

BatchDRK-TES-JUL26

Sample description >

See picture or pictures attached.

Tests requested >

Tesamorelin analysis

Results >



Tesamorelin	9.83 mg
Purity	97.259%

Comments >

Analysis conducted > 29 JUL 2026

Signature > 

Verify this test at www.janoshik.com/verify/ with the following unique key

14RSNWH7I5DF




SHIMADZU
LabSolutions


Analysis Report

<Sample Information>

Sample Name	: 26640-1 Tesa 10mg +LCMS scr	Sample Type	: Unknown
Sample ID	: 124		
Data Filename	: 26640-1 Tesa 10mg +LCMS scr_084.lcd		
Method Filename	: Peptidy_korekce_SIF.lcm		
Batch Filename	: 28_07_2026_01.lcb		
Vial #	: 4-30		
Injection Volume	: 1 uL		
Date Acquired	: 29.07.2026 4:32:49	Acquired by	: System Administrator
Date Processed	: 29.07.2026 12:46:59	Processed by	: System Administrator

<Chromatogram>

mAU



<Peak Table>

PDA Ch2 280nm

Peak#	Ret. Time	Area	Conc.	Unit	Name	Area%
1	3.689	2182	0.00			2.216
2	3.742	95741	9.83	mg	Tesamorelin	97.259
3	4.022	517	0.00			0.525
Total		98439				100.000

Task Number #210098

Testing ordered > 19 JUL '26

Sample received > 23 JUL '26

Clientdrakonpeptides.com

SampleTesamorelin 10mg

ManufacturerUnknown

BatchDRK-TES-JUL26

Sample description >

See picture or pictures attached.

Tests requested >

Endotoxin Analysis

Results >



Endotoxin	1.831 EU/Vial.

Comments >

Analysis conducted > 30 JUL 2026

Signature > 

Verify this test at www.janoshik.com/verify/ with the following unique key

HZBHWR2FM2E3

Task Number #210099

Testing ordered > 19 JUL '26

Sample received > 23 JUL '26

Client drakonpeptides.com

Sample Tesamorelin 10mg

Manufacturer Unknown

Batch DRK-TES-JUL26

Sample description >

See picture or pictures attached.

Tests requested >

LCMS Screening For Peptide Contamination

Results >



LCMS screening	

Comments >

LCMS analysis identified the presence of mannitol, tesamorelin, and degradation products. LCMS spectra are attached. No significant contamination detected.

Analysis conducted > 28 JUL 2026

Signature > 

Verify this test at www.janoshik.com/verify/ with the following unique key

HP2UGY9E9TPK




JANOSHIK

==== Shimadzu LabSolutions Data Report ====

<Spectrum>

Line#:1 R.Time:0.254(Scan#:119)
MassPeaks:30
RawMode:Single 0.254(119) BasePeak:205(2371532)
BG Mode:None Segment 1 - Event 1




**JANOSHIK**

==== Shimadzu LabSolutions Data Report ====

<Spectrum>

Line#:1 R.Time:1.486(Scan#:692)
MassPeaks:59
RawMode:Single 1.486(692) BasePeak:857(970064)
BG Mode:None Segment 1 - Event 1




**JANOSHIK**

==== Shimadzu LabSolutions Data Report ====

<Spectrum>

Line#:1 R.Time:1.455(Scan#:678)
MassPeaks:37
RawMode:Single 1.455(678) BasePeak:734(677636)
BG Mode:None Segment 1 - Event 1





JANOSHIK

==== Shimadzu LabSolutions Data Report ====

<Chromatogram>

