ANALIZA BIAŁEK

Scan for verify

100019227

Laboratory of Analysis

Client: Drakon Peptides
Order number: 100019227
Name: DRAKON PEPTIDES RETATRUTIDE
Number of samples: 1 vial

Batch number: DRK-R40-OCT26
Receive date: 2026-10-07
Start date: 2026-10-12
End date: 2026-10-26

Sample identification

Substance	Method of taking a sample	Form
Retatrutide	Sample provided by customer	White powder in a glass vial

Results

Substance	Sample name	Content	Purity
Retatrutide	DRAKON PEPTIDES RETATRUTIDE 40MG	41,08mg	99%
Molecular Mass	DRAKON PEPTIDES RETATRUTIDE 40MG	4731,25m/z	------
Endotoxin Level	DRAKON PEPTIDES RETATRUTIDE 40MG	7,87EU/vial	------

Research comment

The test sample was applied in the amount of 3µl to the BioBasic 8 Thermo Scientific column. During the course of the analysis with the method a linear gradient in the solutions of H2O and Acetonitrile with the addition of 0.1% trifluoroacetic acid, a peak of the sought substance appeared. Its content, when compared to the reference chromatogram, is 41,08mg. Height the measurement error is +-0,05mg. The actual state of the sample sent before the test is presented on the photo below.

Photo of the test sample



Password

YGK38TIF

Analiza Bialek sp. z o.o.
KRS: 0000966168
NIP: 8992922970
REGON: 521751875

Headquarters and Laboratory:
ul. Duńska 9, 54-427 Wrocław
Budynek DELTA


Lab technician:
Analysis and research specialist
mgr Arkadiusz Zając

ANALIZA BIAŁEK

Scan for verify

100019227

Chromatogram



Method

Time	H2O + 0,1% TFA	ACN + 0,1% TFA
0min	95%	5%
2min	95%	5%
12min	5%	95%
15min	5%	95%

Password

YGK38TIF

Analiza Bialek sp. z o.o.
KRS: 0000966168
NIP: 8992922970
REGON: 521751875

Headquarters and Laboratory:
ul. Duńska 9, 54-427 Wrocław
Budynek DELTA

Lab technician:
Analysis and research specialist
mgr Arkadiusz Zając

ANALIZA BIAŁEK

Scan for verify

100019227

Laboratory of Analysis

Client: Drakon Peptides
Order number: 100019227
Name: DRAKON PEPTIDES RETATRUTIDE
Number of samples: 1 vial

Batch number: DRK-R40-OCT26
Receive date: 2026-10-07
Start date: 2026-10-12
End date: 2026-10-26

Sample identification

Substance	Method of taking a sample	Form
Retatrutide	Sample provided by customer	White powder in a glass vial

Results

Test	Sample name	Result
TYMC	DRAKON PEPTIDES RETATRUTIDE 40MG	PASS
TAMC	DRAKON PEPTIDES RETATRUTIDE 40MG	PASS
Presence of presumptive Escherichia Coli	DRAKON PEPTIDES RETATRUTIDE 40MG	NOT DETECTED
Presence of Salmonella Spp.	DRAKON PEPTIDES RETATRUTIDE 40MG	NOT DETECTED

Heavy metal	Sample name	Result
Arsenic	DRAKON PEPTIDES RETATRUTIDE 40MG	<0,010mg/kg
Cadmium	DRAKON PEPTIDES RETATRUTIDE 40MG	0,103mg/kg
Lead	DRAKON PEPTIDES RETATRUTIDE 40MG	0,335mg/kg
Mercury	DRAKON PEPTIDES RETATRUTIDE 40MG	<0,010mg/kg

Password

K0JZHAHF

Analiza Bialek sp. z o.o.
KRS: 0000966168
NIP: 8992922970
REGON: 521751875

Headquarters and Laboratory:
ul. Duńska 9, 54-427 Wrocław
Budynek DELTA

Lab technician:
Analysis and research specialist
mgr Arkadiusz Zając

ANALIZA BIAŁEK

Scan for verify

100019227

Photo of the test sample



Password

K0JZHAHF

Analiza Bialek sp. z o.o.
KRS: 0000966168
NIP: 8992922970
REGON: 521751875

Headquarters and Laboratory:
ul. Duńska 9, 54-427 Wrocław
Budynek DELTA

Lab technician:
Analysis and research specialist
mgr Arkadiusz Zając